HYPERSPHERE

BY HYPERICE

OPERATING INSTRUCTIONS | PLEASE READ PRIOR TO USE

The HYPERSPHERE is a cutting-edge fitness and recovery device that uses pressure and vibration to improve the body's overall performance. The HYPERSPHERE is manufactured using industry best components providing just the right amount of pressure and vibration. The HYPERSPHERE offers people of all ages and levels the combined benefits of compression (pressure) and vibration to improve the body's overall performance, enhance recovery, flexibility, range of motion to prevent injury and thus improve overall well-being.

Charging:

- To turn on, press down on the power button for three seconds until LED lights illuminate. Press again once for level one, twice for level two and three times for level three. The blue LED lights will indicate the corresponding power level.
- Do not use directly on a hard surface. Use only on a carpeted floor or a fitness mat.

Operating:

- With Power Switch in ON position, press the power button once to turn unit on to level one, twice to turn unit to level two and three times for level three. Pressing the power button again will turn the unit off. The blue LED lights will indicate the corresponding power level.
- The unit will turn off after 10 minutes of use. To continue use, simply press the power button again to the desired level.
- To turn off, press down on the power button for three seconds. The unit is fully turned off when all the LED lights are off.
- The unit will turn off after 10 minutes of use. To continue use, turn unit on again and set to desired level.

Maintenance , Cleaning and Storage:

- Use a towel or soft cloth to wipe any moisture from unit.
- Keep unit away from heat sources, liquids, and children.

Warnings and Caution

- Do not immerse unit in water. Keep liquids away from ventilation ports, buttons and charging port.
- Do not remove screws or attempt to disassemble.
- Do not operate continuously for more than one hour. After one hour use, allow unit to rest for 30 minutes before re-using.
- Unplug the unit after charging or prior to use.
- Do not use directly on a hard surface. Use only on a carpeted floor or a fitness mat.
- Persons weighing above 250 lbs /113 kgs may damage or break the unit during normal intended use.
- This is not a toy. For adult use only. Do not use if injured. Consult your doctor before using this product.
- Use product only according to Instructional Use materials available at www.hypence.com.


DISCLAIMER

Please do not use the HYPERSPHERE, or any vibration high intensity exercise device without first obtaining approval from your doctor if any of the following apply:

Pregnancy, diabetes with complications such as neuropathy or retinal damage, wear of pace-makers, recent surgery, epilepsy or migraines, herniated disks, spondylolisthesis, spondylolysis, or spondylosis, recent joint replacements or IUD's, metal pins or plates or any concerns about your physical health. Frail individuals and children should be accompanied by an adult when using any vibration device.

These contra indications do not mean that you are not able to use a vibration or exercise device but we do advise you to consult a doctor first.

Ongoing research is being done to the effects of vibration exercise for specific medical disorders. This will likely result in the shortening of the list of contra indications as shown above. Practical experience has shown that integration of vibration exercise into a treatment plan is advisory in a number of cases. This must be done on the advise of and in accompaniment with a doctor, specialist or physiotherapist.

1. Charging Port
2. Power Button
3. Battery level indicator LED (green, red)
4. Frequency level indicator LED (blue)
5. AC wall charger (16VDC-1A)

Specifications

- AC Charger: 110-240V 50/60Hz - .5A - 16VDC - 1A
- Rechargeable Lithium Ion battery 2600mah
- Frequency: Low 45Hz, Medium 68Hz, High 92Hz
- Weight: 2.8 lbs / 1.3 kg (including charger)
- Designed in USA. Assembled in China

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